



DOTS Platform

BUILD A SMARTER SHAKE FLASK, ONE PARAMETER AT A TIME

DOTS brings advanced sensing and control options, simple handling, and lower costs than similar fermentation platforms.

The DOTS Platform: The New Status Quo for Shake Flask Fermentations



The DOTS Platform turns standard shake flasks into smart mini bioreactors.

The DOTS Platform: The New Status Quo for Shake Flask Fermentations

Liquid Injection System (LIS)

Automated feeding of liquids into shake flask cultures. LIS enables bioreactor-like options in shake flasks: Fed-batch, DO-based feeding, biomass-based feeding, automated promotor induction, and more.

Multiparameter Sensor (MPS)


Optical and non-invasive monitoring of various process parameters like biomass, fluorescence, dissolved oxygen, and many more.

DOTS Software

Powerful DOTS Software for easy sensor handling and real-time data visualization.

Dissolved Oxygen (DO) Sensor Pill

Each pill contains a DO-sensitive chemosensor that is read out by the Multiparameter Sensor (MPS).



Choose from an array of parameter-specific sensor-and-software modules to build the smart shake flask that meets your needs.


The DOTS Platform: The New Status Quo for Shake Flask Fermentations



Biomass



Dissolved Oxygen (DO)



Fluorescence



**Parameter-based
Feeding**



**Environmental
Parameters**

With the DOTS Platform, you can continue to expand your shake flask sensing capabilities as new sensors are added in the future.



Sensing & Controlling

DOTS combines sensors and actuators with a powerful software to create bioreactor-like shake flasks with actionable insights into your bioprocess.

DOTS Platform Hardware & Software Components

Multiparameter Sensor (MPS), Adapter & USB Hub



The MPS is an optical sensor, capable of reading out various signals from one shake flask culture. It is positioned in the adapter, underneath the shake flask, allowing for increased stability. The USB Hub bundles data from multiple MPS and forwards it to the DOTS Software.


Parameter Specific Sensors & Actuators



Combine with Dissolved Oxygen (DO) Sensor Pills for online dissolved oxygen (DO) monitoring in shake flasks.

The Liquid Injection System (LIS) enables the automated feeding of liquids into shake flask cultures.

DOTS Software



DOTS Software enables a simplified control of sensors and visualizes the received data from all monitored shake flasks in real-time.

DOTS current sensors and actuators are designed to simplify bioprocessing for a variety of applications.

The Current DOTS Portfolio



**Multiparameter
Sensor (MPS)**



**Dissolved Oxygen
(DO) Sensor Pills**



**Liquid Injection
System (LIS)**



The MPS is an optical sensor technology for the monitoring of multiple parameters in shake flasks.

Multiparameter Sensor (MPS)

Check Out
Individual
Product Presentation



a platform by S&B



Key Facts

Most versatile shake flask sensor on the market

Turn your shake flask into a low-cost, high-throughput mini bioreactor

Monitor a broad range of Critical Process Parameters (CPPs)

Biomass, Dissolved Oxygen (DO), Fluorescence, and more!

Save time with automated, online monitoring

Installed underneath the shake flask and measures through the vessel wall

Works on a wide range of shake flasks

Compatible with shake flasks from 100 mL - 2000 mL, clamps and Sticky Mats

For a broad range of organisms

Bacteria, (filamentous) fungi, algae, archaea, and plant cells

Powerful DOTS Software

Simplified sensor control and data visualization for improved comparability

First pill-based optical sensor for online dissolved oxygen (DO) monitoring in shake flasks.

Dissolved Oxygen (DO) Sensor Pills

Check Out
Individual
Product Presentation



DOTS
a platform by SBI



Key Facts

Novel, patented pill technology

Enables automated, online DO monitoring in shake flasks

Single-use pill design

Factory-calibrated and pre-sterilized for immediate use

Unique pill identification algorithm

Removes the need for sensor alignment

Easy to use

Drop & Go: Easy handling and fast experiment setup

Enhanced control options

Combine with LIS for DO-based feeding in shake flasks

Powerful DOTS Software

Simplified sensor control and data visualization for improved comparability

The LIS is the first technology allowing for automated liquid feeding in shake flasks.

Liquid Injection System (LIS)

Check Out
Individual
Product Presentation
 DOTS
a platform by sbi



Key Facts

Wireless control and monitoring

Control and monitor your feeding experiments wirelessly with the DOTS Software

New application: parameter-based feeding

Combine with MPS and enable biomass-based or DO-based feeding in shake flasks

Compatible with various substances

Sugars (e.g., 40% glucose), alcohols (e.g., 50% methanol), glycerol, acids, bases, ...

Create any feeding profile you want

Single or multi shot, exponential, or constant feeding

Easy to install and use

Fill the sterile cartridge, program the LIS drive and start feeding your culture

Powerful DOTS Software

Simplified sensor control and data visualization for improved comparability



Sensor Handling & Data Visualization

The DOTS Software enables easy sensor handling and experiment set up.

Exemplary Screenshots

Create an experiment with pre-defined application templates

The screenshot displays the DOTS software interface for creating an experiment. The top navigation bar includes the DOTS logo, a search bar, and a user profile icon. The main content area is divided into two sections: 'Basic Settings' (Step 1) and 'Device Assignment' (Step 2).

Basic Settings:

- Basic information:** Includes radio buttons for 'Quick start application templates' (selected) and 'Custom application templates'. A dropdown menu shows 'Biomass monitoring <150...'. There is also a 'Select template' dropdown.
- Enable template configuration step
- Experiment name *:** Text input field containing 'Strain B, 2% glucose'.
- Number of objects *:** Text input field containing '1'.
- Project *:** Dropdown menu showing 'Prosugar'.
- Advanced settings

Device Assignment:

- Strain B, 2% glucose:** A central box containing a tree structure of tasks:
 - Biomass monitoring <150 rpm (Shake flask):** A blue box at the top.
 - Cell growth (CGQ):** A box with a 'Planned' status.
 - OD600 (Offline):** A box with a 'Planned' status.

Assign sensors to planned experiments via drag and drop

The screenshot displays the DOTS software interface for assigning sensors to a planned experiment. The top navigation bar includes the DOTS logo, a search bar, and a user profile icon. The main content area is divided into four steps: 'Basic Settings', 'General Task Configuration', 'Replicate configuration', and 'Device Assignment' (Step 4).


Device Assignment:

- Experiment:** MET25 induction strain A (0/1).
- Tasks:** A tree structure of tasks:
 - Biomass-based feeding (Shake flask):** A blue box at the top.
 - Feeding (LIS):** A box with a 'Planned' status and a 'No device connected' warning.
 - Biomass monitoring (CGQ):** A box with a 'Planned' status and a 'No device connected' warning.
 - OD600 (Offline):** A box with a 'Planned' status.
- Device List:** A list of devices with status indicators:
 - LIS-0025478:** LIS-LIS0-254869 / 25. Status: 80% (green), Not connected (red), Unassigned (blue).
 - CGQ-SP-02548:** CGQ-8-0025 / Port1. Status: Not connected (red), Unassigned (blue), Free (green).
- Instructions:** 'Drag and drop device in order to connect it to process / task.'
- Buttons:** 'Exit wizard', 'Back', and 'Create'.


The DOTS Software provides a comprehensive overview of your experiments and visualizes your data in real-time.

DOTS Software Modules for Data Visualization

The dashboard schema provides an overview over all running, planned, or finished experiments



Data visualization tools enable a customized display of data in real-time






Applications


Exemplary data: By measuring fluorescence emission of different wavelengths, photoswitch processes can be tightly monitored.

Fluorescence Monitoring: *E.coli* Expression of a Photoswitch Protein



With the DOTS Platform, parameter-based feeding in shake flasks is now a reality.


Parameter-based Feeding With The DOTS Platform



Exemplary data: Feeding with the Liquid Injection System starts when a preset DO-threshold, measured with DO Sensor Pills, is reached.


DO-based Methanol Feed to *Pichia pastoris* Cultures to Keep Promotor Activity Constant

- With a DOTS-integrated controller, methanol feed is adjusted
- As soon as the DO of the culture has recovered (>80% DO) the Liquid Injection System (LIS) starts feeding
- This enables ideal cell viability while keeping promotor activity constant



Exemplary data: With biomass-based feeding, the Liquid Injection System starts at a preset biomass or growth rate threshold.

Biomass-based Feeding: Methanol Induction With *Pichia pastoris*



sbi SCIENTIFIC
BIOPROCESSING

Let's Connect!

insights@scientificbio.com
www.scientificbio.com

@scientific bioprocessing 

@scientific bioprocessing 

@scientific_bio 

