


A NEW APPROACH

A NEW APPROACH TO CELL CULTIVATION: OPTICAL SENSORS


- The most commonly-used sensors in bioreactors are electrochemical probes but are unsuitable for T-flasks and microfluidics devices.
- Therefore researchers suffer from a lack of fundamental knowledge of the physical background and controlling parameters in early stage bioprocessing development.
- There is now an acute need for improvements. Real-time on-line monitoring of key parameters in bioprocessing scouting devices is needed to produce conditions representative of manufacturing.¹

HYBRIDOMA CELL GROWTH STUDY IN ROCKING VS STATIC T-FLASKS



- Non-adherent SP2/O-based mouse hybridoma cells secreting IgG3
- T-75 flasks, static and on a rocker platform

Cell Culture Comparability Study Between Static and Rocking T75-flasks After Cell Vial Thaw



- Cell cultures grew faster and more densely in rocking than static t-flasks
- Rocking flasks yielded higher viable cell density (VCD) (1.8-fold increase) and 25% less lactate production compared to static flasks.
- DO in the rocking flask never dropped below 70% versus the static flask dropping to zero in two days leaving them anoxic for hours at a time.
- Agitation had a positive impact on maintaining physiological pH
- 31% higher antibody titer than static flasks
- Static T-flasks should be passaged by day 2 where rocking flasks never reached 0% DO and can be passaged by 3.5 days with cell viability profile remaining above 90%.

- Hypoxia can drive and maintain genetic instability, and hypoxic cells can acquire a mutator phenotype that consists of decreased DNA repair, an increased mutation rate and increased chromosomal instability.²
- Monitoring DO profile data may increase cell line genetic stability.³


ID-SENSORS

SCALE DOWN STUDY

SBI's low cost, disposable DO and pH sensors in rocking T-flasks matching $k_L a$ with a 10L wave bioreactor¹


- Compared three T-flask geometries (T-25, T-75 and T-150) with optical ID-Sensors pH affixed to the flasks.
- Measured $k_L a$ of spinner flasks and found the T-flask superior.
- The ID-Sensor DO was used to accurately estimate the $k_L a$ in rocking T-flasks allowing researchers to match it to a Cultibag wave bioreactor.
- T-flasks have high surface area-to-volume ratio (A/V , cm^2) one order of magnitude higher than other process scouting devices (PSDs) even at low Reynolds values. $k_L a$ increases as A/V increases.¹
- Traditionally, this was thought to compensate for lack of agitation.
- Traditional wisdom has been shown to be incorrect.

$k_L a$ Comparability Between Rocking T25-flasks and Other Devices



- Rocking T-Flasks Have Superior Oxygen Transfer Capabilities Versus Other PSDs And Higher $k_L a$ At Lower Re Values.

OPTICAL SENSING



- Fluorescence-based technology
- Intelligent, dynamic (ID) products
- ID-Reader sits outside of culture vessel
- ID-Sensor is affixed inside of culture vessel
- Vessel wall must be transparent and transmit light $>320\text{nm}$

CONCLUSION

Four features required to deliver the power to optimize cell culture:



- Reliable sensors for real-time monitoring from smallest to largest vessel
- Automatic with control and closed loop feedback
- Relevant physiological parameters at the cellular level
- Economically remarkable and easy to demonstrate ROI

CONTACT US

Jake Boy, Senior Applications Scientist
jboy@scientificbio.com

REAL TIME MONITORING

ID-READER AND ID-FIBER OPTIC READER



ID-DEVELOPER'S KIT


- Real-time, non-invasive monitoring in one kit
- Optical sensors for pH and dissolved oxygen (DO), readers and software



ID-ROCKER

Rocking and Instrumented T-Flasks are Better Flasks

- Adjustable T-flask rocking speed
- Closed loop control of rocking speed with feedback from the sensor measurements



ID-DATA HUB

Software for Optical Sensors

- Select the type of sensor you are using, pH or dissolved oxygen
- Select data collection mode
- Watch live data display
- Record data
- Single mode to capture single data points
- Scan mode to gather data at 10-second or higher intervals



BIBLIOGRAPHY

- Vallejos JR, Micheletti M, Brorson KA, Moreira AR, Rao G. Optical Sensor Enabled Rocking T-Flasks as Novel Upstream Bioprocessing Tools. *Biotechnology and Bioengineering*. 109(9):2295 2012
- Bristow RG, Hill RP. Hypoxia, DNA repair and genetic instability. *Nature Reviews Cancer* 8(3):180 2008
- Vallejos JR, Brorson KA, Moreira AR, Rao G. Dissolved Oxygen and pH Profile Evolution After Cryovial Thaw and Repeated Cell Passaging in a T-75 Flask. *Biotechnology and Bioengineering* 105(6):1040 2009