

CELL GROWTH QUANTIFIER BioR (CGQ BioR)

Biomass Monitoring For Bioreactors





Cell Growth Quantifier BioR (CGQ BioR)

Biomass Monitoring For Bioreactors

Key Facts

✓ **For a broad variety of organisms**

Bacteria, (filamentous) fungi, algae, archaea, and plant cells

✓ **Compatible with various scales of bioreactors**

Mini bioreactor systems, benchtop bioreactors, production scale, and more

✓ **A cost-effective solution**

An affordable technology for one of the key bioprocess parameters

✓ **Flexible integration**

Integration of biomass data into other software available via analog, OPC, and other interfaces

✓ **Compatible with various scales of bioreactors**

Install the sensor to the outside of the bioreactor, without blocking bioreactor ports

✓ **DOTS Software**

Powerful software for easy sensor handling and real-time data visualization